

## **Third class (PCR and digestion, 2006)**

### **PLAN FOR THE DAY**

- a) Set up PCR reaction of GFP and potassium channel
- b) Check the size of expected PCR product by computer
- c) Check PCR product on gel and take picture
- d) Ethanol precipitate DNA
- e) Set up digestion of PCR product and vector.

### **1. PCR reaction**

Because PCR takes some time, we will start the reaction first. As you come in, please set up PCR reaction.

Mix the following solution in a thin-wall PCR tube (200 µl).

- 82 µl of water
- 10 µl of buffer [for taq Polymerase]
- 2 µl of CaMKII 5' primer (50 pmol/µl)
- 2 µl of CaMKII 3' primer (50 pmol/µl)
- 2 µl of dNTP (10 mM each)
- 1 µl of Taq polymerase
- 1 µl of CaMKII plasmid DNA (100 ng/µl)

Mix by gently tapping.

Spin for 3 seconds in small centrifuge

Split the contents of the tube into 3 thin-walled tubes with 33µl each.

Start PCR reaction at the following condition (the machine will be preset to this)

1. 94 °C, 1 min
  2. 55 °C, 0.5 min
  3. 72 °C, 1 min
- Repeat the same reaction for 20 times

### **2. Estimate the size of PCR product**

The primer sequences are as follows;

```
5' GAT CCT CGA GCT ATG GCT ACC ATC ACC TGC ACC
3' GAT CG AAT TCG CAT GCC CTG GCC GTT GCC TTC AAT G
```

We are amplifying rat CaMKII alpha coding region..

Using these hints and the DNA database, estimate the size of expected PCR product.

### **3. Run PCR product on agarose gel**

(prepare an agarose gel as soon as you are done with the PCR)

Take 5  $\mu$ l from the tube, mix with 1  $\mu$ l of DNA loading dye (Xc for xylene cyanol) on parafilm, and load on to a gel.

Load 4  $\mu$ l of 1 kb DNA ladder

Once everything is loaded, set voltage to 150 volts and start running the gel

Run gel for 15 minutes

Stop gel and take image with UV light

#### **4. Ethanol precipitate PCR product**

You can do this while running the gel.

To the rest (95  $\mu$ l) of PCR product, add the following

10  $\mu$ l of sodium acetate 3 M

250  $\mu$ l of 100% ethanol

Mix well by shaking

Spin in big centrifuge at maximum speed (14000 rpm) for 5 minutes

Discard supernatant in waste bucket

Add 1 ml of 70% ethanol

Discard supernatant in waste bucket

Air dry tube for 5 minutes

Resuspend in 54  $\mu$ l of water

#### **5. Set up the digestion of vector and insert**

For the vector: (Note: make sure that you take 1  $\mu$ l or less of enzyme!)

6  $\mu$ l of pEGFP-C1

6  $\mu$ l of EcoRI buffer

46  $\mu$ l of water

3  $\mu$ l of BSA

3  $\mu$ l of EcoRI

3  $\mu$ l of XhoI

1  $\mu$ l of shrimp alkaline phosphatase

Finger-mix, spin for 5 seconds in small centrifuge.

For the PCR insert:

Take 54  $\mu$ l of the PCR product

Add 6  $\mu$ l of EcoRI buffer

3  $\mu$ l of EcoRI

3  $\mu$ l of XhoI

Finger-mix, spin for 5 seconds.

Incubate both reactions with insert and vector at 37C for at least three hours.

We will continue next week.